


# Cyclic Compression Porometer


## Principles of Operation

The instrument applies cyclic stress on a sample and measures the pore structure characteristics after a desired number of cycles. The sample is loaded in the sample chamber and is subjected to stress cycles in the specified stress limits. At the end of the desired number of stress cycles, the pressure of a nonreacting gas on one side of the sample is increased to initiate gas flow through pores. The gas pressure and flow rates are measured. The pores in the sample are spontaneously filled by a wetting liquid. The gas pressure and flow rates are measured through the wet sample. After acquisition of data, the sample is re-wetted and again subjected to cyclic compression. Pressure and flow rates are measured after the desired number of cycles. The test is continued to acquire data as a function of number of stress cycle.



## Features

- Fully Automated
- After desired number of stress cycles
  - automatically interrupts analysis
  - performs tests
  - acquires data
  - continues to cyclically stress the sample
- Very little operator involvement
- Operator adjustable
  - stress limits in a cycle
  - number of cycles after which data is to be acquired
  - stress free duration for data acquisition
  - total number of cycles
- Concurrent measurement of compressive strain in the sample as a function of stress cycle
- Windows based simple operation



## Specifications

**Pore Size Range:**  
0.013 - 500 microns

**Permeability Range:**  
1 x 10<sup>-8</sup> - 50 darcies

**Sample Size:**  
1.75" - 2.5" diameter

**Pressure Range:**  
0 - 500 psi

**Pressurizing Gas:**  
Clean, dry, and compressed air or  
nonflammable and noncorrosive gas

**Pressure Transducer Range:**  
0 - 500 psi

**Resolution:**  
1 in 20,000

**Accuracy:**  
0.15% of reading

**Mass Flow Transducer Range:**  
10 cm<sup>3</sup>/min - 500,000 cm<sup>3</sup>/min

**Power Requirements:**  
110/120 VAC, 50/60 Hz  
(Others Available)

**Dimensions:**  
30" H x 19" W x 18.5" D

**Weight:**  
100lbs

## Effects of Cyclic Compression on Pore Diameter of Felts

Material	Maximum Compressive Stress, psi	# of Cycles	% Change in Bubble Point	% Change in Mean Flow Pore Diameter
Felt #1	500	15	-71.1	-30.3
Felt #2	750	2000	-88.4	-15.8

## Other Products

Average Fiber Diameter Analyzer  
Bubble Point Tester  
Capillary Flow Porometer  
Capillary Condensation Flow Porometer  
Complete Filter Cartridge Analyzer  
Clamp-On Porometer  
Compression Porometer  
Custom Porometer  
Cyclic Compression Porometer  
Envelope Surface Area Analyzer  
Filtration Media Analyzer  
High Flow Porometer  
Integrity Analyzer

In-Plane Porometer  
Microflow Porometer  
Nanopore Flow Porometer  
QC Porometer  
Diffusion Permeameter  
Gas Permeameter  
Liquid Permeameter  
Vapor Permeameter  
Water Vapor Transmission Analyzer  
Liquid Extrusion Porosimeter  
Mercury/Nonmercury Intrusion Porosimeter  
Vacuapore  
Water Intrusion Porosimeter (Aquapore)

BET Liquisorb  
BET Sorptometer  
Gas Pycnometer  
Mercury Pycnometer

**Also Available:**  
Testing Services  
Consulting Services  
Short Courses

**Buy Rent Lease**